

Ramister-C

Isterometro curvo
Curved Hysterometer



- **Ramister-C:** Fornisce una misura approssimativa della profondità della cavità uterina. Ricurvo e graduato con scala da 5 a 11 cm (risoluzione 0,5 cm)* con lettura diretta sull'impugnatura. Facile posizionamento grazie al fermo scorrevole da collocarsi sull'orifizio cervicale. Leggermente conico di diametro da 3 a 4 mm.

* Fornisce una misura approssimativa.

- **Ramister-C:** It provides an approximate measure of the depth of the uterus. Curved and graduated, ranging from 5 to 11 cm (0.5 cm resolution)* with direct reading on the handle. It is easy to place, thanks to the adjustable positioner to be set on the cervix external orifice. Slightly conical configuration (3 to 4 mm o.d.).

* It gives an approximate measurement.

Cod.	Descrizione / Description	Pz/Scatola - Pcs/Box
720001	Lunghezza 24,5 cm Length 24.5 cm	50 pz./pcs.

Fornito in confezione singola - monouso - sterilizzazione a EtO - marcato CE 0051

Supplied in single packaging - disposable - EtO sterilization - CE 0051 marked

