




**THE BEST QUALITY
IN THE WORLD,
TO THE WORLD.**



**iTiPremium
HAND INSTRUMENTS**

Products Name		Size	Taper	Working Length	Material	Pack
NiTi FLARE FILES		#15-40	0.05	21mm, 25mm	NiTi	6pcs
NiTi H-FILES		#15-40	0.02	21mm, 25mm	NiTi	6pcs
NiTi FLARE SPREADERS		#15-40	0.05	21mm, 25mm	NiTi	6pcs
NiTi PLUGGERS		#15-40	0.035	21mm, 25mm	NiTi	6pcs

8-3 KIOHARA INDUSTRIAL PARK, UTSUNOMIYA TOCHIGI, 321-3231 JAPAN

Phone: +81-28-667-8592 FAX: +81-28-667-8593

E-mail: dental.exp@ms.mani.co.jp

Website: <http://www.mani.co.jp/en/index.html>

MANI®