



## MAIN APPLICATIONS

Emergency Dept.

Home-Care

Patient Transfer

Room / Ward

Tracheostomy

Veterinary

Internal Rechargeable Battery

ASKIR 36BR is a portable suction unit for professional use and with a powerful aspiration of max 36 l/min, conceived to be lightweight for easy transport. Three different options for operation: AC/DC adapter, internal lead or lithium-ion rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level.

The newly-designed ergonomic handle ensures a comfortable and secure grip.

## ACCESSORIES INCLUDED

- Liquid Collection Jar with overflow valve system (different options, see below)
- Silicone Tubes  $\varnothing$  8x14mm (autoclavable) - Patient Tube length 150cm
- Conical Connector  $\varnothing$ 10-11-12mm
- Antibacterial & Hydrophobic Filter (single-patient)
- Internal rechargeable Lead Battery
- Universal switching power adapter and Cable + 12V Car Adapter

## OPTIONS

- Stand on 5 Castors  $\varnothing$  50mm with Brakes (antistatic upon request) REF 27730
- Carrying Bag (for 1000ml jar) SP 0207/02
- Carrying Bag (for 2000ml jar) SP 0207
- Wall Mount Brackets 27736
- Vacuum regulator footswitch /

## TECHNICAL FEATURES

Motor	Oiless and maintenance-free piston pump
Power Feeding	14V $\approx$ 4A by AC/DC Universal Adapter 100-240V / 50-60Hz - 100VA (model UE60-140429SPA1)
	12V $\approx$ 4A Lead rechargeable battery - Autonomy 60 min. - Recharging Time 4 hours
	12V $\approx$ 4A Car adapter
ON/OFF Switch	Mechanical switch
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
IP Code	IP21
Max Vacuum (adjustable)	-0.80 bar -80 kPa -600 mmHg
Max free air flow rate	36 l/min
Noise Level	65,5 dBA (with jar) / 68,5 dBA (without jar)
Duty cycle	Non-stop operation (by mains only)
Weight	kg 5.31 (unit packed with all accessories)
Size	cm 35 x 18 x 21
Years of Warranty	2
Shipping carton	3
Place of Manufacturing	Italy

Italian  
Medical  
Touch  
MADE IN ITALY

## AVAILABLE CONFIGURATIONS

RE 410200/03	ASKIR 36BR with 1000 ml autoclavable jar in polycarbonate Makrolon 2858
RE 410200/09	ASKIR 36BR with 2000 ml autoclavable jar in polycarbonate Makrolon 2858
RE 410201	ASKIR 36BR with 1000 ml Flovac® disposable liner and reusable container
RE 410201/01	ASKIR 36BR with 2000 ml Flovac® disposable liner and reusable container
	Jars in APEC® autoclavable up to 143°C available upon request
	Disposable jar 1000 ml available upon request

CE 0123