

## ST-C/ST-D Mechanical Operating Table



ST-C

Technical data	
Dimensions of table top(mm)	1960*480
Upward folding angle of head section	35°
Downward folding angle of head section	40°
Downward folding angle of leg section	90°
Upward folding angle of back rest	75°
Reversed Trendelenburg of table	20°
Trendelenburg of table	15°
Left inclining angle	15°
Right inclining angle	15°
Lift of waist section(mm)	80
Minimum height of table(mm)	720
Lift of table(mm)	200

Standard Accessories	
Mattress	1 set
Leg rest	2 pcs
Arm rest	2 pcs
Shoulder support	2 pcs
Anesthetic screen	1 set
Body restraining strap	1 pc

ST-C&ST-D mechanical operating tables can meet different surgical requirements and can be used for general and specialized surgery.  
ST-D Mechanical Operating Table is the same as ST-C except the difference of leg section.

- ST-C: One piece leg section.
- ST-D: Divided leg section

### Optional Accessory

Waist bridge	1set
--------------	------



ST-D

## ST-3008/ST-3001 Series Mechanical Operating Table

ST-3008 series operating table is necessary for various operation which is designed for operation on head neck abdomen chest the four limbs and others in general surgery. Using the current latest oil pump for hydraulic elevation, and the other operations are controlled manually

- ST-3008A The table base and pump shell with stainless steel
- ST-3008B The table base, top surface and pump shell with stainless steel
- ST-3008C The table base, pump shell with stainless steel, two leg plates, double-layer desk for X-RAY
- ST-3008D The table base, pump shell with stainless steel two leg plates



ST-3008D



ST-3008C

### ST-3008 Series

Technical data	
Table height (mm)	800~1040
Length and width (mm)	2100*480
Reversed trendelenburg	35°
Trendelenburg	20°
Left inclining angle	20°
Right inclining angle	20°
Upward folding angle of back rest	75°
Downward folding angle of back rest	20°
Upward folding angle of head rest	90°
Downward folding angle of head rest	90°
Downward folding angle of leg rest	90°

ST-3001 series operating table is necessary for various operation which is designed for operation on head neck abdomen chest the four limbs and others in general surgery. Using the current latest oil pump for hydraulic elevation, and the other operations are controlled manually

- ST-3001A The table base and pump shell with stainless steel
- ST-3001B The table base, top surface and pump shell with stainless steel

### ST-3001 Series

Technical data	
Table height (mm)	800~1040
Length and width (mm)	2100*480
Reversed trendelenburg	35°
Trendelenburg	20°
Left inclining angle	20°
Right inclining angle	20°
Upward folding angle of backrest	75°
Downward folding angle of backrest	20°
Upward folding angle of headrest	90°
Downward folding angle of headrest	90°
Downward folding angle of legrest	90°

### ST-3008 / ST-3001

Standard Accessories	
Anesthetic screen	1pc
Leg rest	2pcs
Shoulder support	2pcs
Waist rest	2pcs
Arm rest	1pc



ST-3001B