

ST-A Mechanical Operating Table



- The lifting of the table is driven by oil pump
- The tenderlenburg, lateral tilt left/right, backplate upwa by air cylinder
- The shell of the base is stainless steel
- The mattress are made of PU material. it is soft and comfortable, with plasticity against pressure ulcers effectively

| Technical Data | |
|-------------------------------------|----------|
| Table height (mm) | 710-910 |
| Length and width (mm) | 1950*480 |
| Reversed trendelenburg | 20° |
| Trendelenburg | 25° |
| Left inclining angle | 15° |
| Right inclining angle | 15° |
| Upward folding angle of back rest | 85° |
| Upward folding angle of head rest | 40° |
| Downward folding angle of head rest | 50° |
| Downward folding angle of leg rest | 90° |
| Lifting height of waist bridge (mm) | 120 |

| Standard Accessories | |
|------------------------|-------|
| Anesthetic screen | 1 pc |
| Leg rest | 2 pcs |
| Body restraining strap | 1 pc |
| Arm rest | 2 pcs |
| Optional Accessory | |
| Memory Foam Mattress | 1 set |



ST-B Mechanical Operating Table



ST-B mechanical Operating Table can meet the requirements for all operating, the operation table can take X ray C-arms image for whole position

- The table top can be rotated 360°
- The shell of the base is stainless steel
- The lifting of the table is driven by oil pump
- The tenderlenburg, lateral tilt left/right, backplate upward/downward, head rest is controlled by air cylinder
- The mattress are made of PU material. it is soft and comfortable, with plasticity against pressure ulcers effectively

| Technical data | |
|-------------------------------------|----------|
| Table height(mm) | 680-880 |
| Length and width (mm) | 1910*500 |
| Reversed Trendelenburg | 20° |
| Trendelenburg | 15° |
| Left inclining angle | 17° |
| Right inclining angle | 17° |
| Upward folding angle of back rest | 75° |
| Upward folding angle of head rest | 35° |
| Downward folding angle of head rest | 40° |
| Downward folding angle of leg rest | 90° |
| Split angle of leg rest | 180° |
| Table-top rotated angle | 360° |
| Lifting height of waist bridge (mm) | 120 |

| Standard Accessories | |
|------------------------|-------|
| Anesthetic screen | 1 pc |
| Leg rest | 2 pcs |
| Body restraining strap | 1 pc |
| Arm rest | 2 pcs |
| Optional Accessory | |
| Memory Foam Mattress | 1 set |