

Cod. Art. :
SL2110



Cradle adjustable height on castors with Trend and reverse-Trend

Height adjustable cot through gas spring control. "X" shaped basement in order to increase the stability and to place the cot as closer as possible to the hospital beds.

Fixed to the central column equipped with twins wheels, two with brake.

Column made of extruded aluminium.

Swinging support system with absolute precision adjustment that allows Trendelenburg of -13° and reverse-Trendelenburg of 13° .

Basin made of non-toxic material, transparent and thermoshaped, provided with place for label and holes at the bottom for air circulation.

Can be provided with box for clothes or personal belongings.

Accessory:



SZ0095-LOWER BASKET FOR SL2110 BASSINET