

TCS-900

OPHTHALMIC UNIT



Features

- Table-top can be opened to facilitate patients' access.
- Two devices can be placed on a slideable table-top
- Motorized up/down chair.
- Phoropter arm can adjust the height and rotation electrically.

Specification

Table	Dimensions	950(W)*853(D)*1908(H)mm
	Table-top size	950(W)*480(D)mm
	Table-top translation range	390mm
	Phoropter arm lifting range	300mm
	Drawer size	535(W)*375(D)*50(H)mm
	Power supply	AC 220V/50Hz, AC 110V/60Hz
	Power consumption	100VA
Chair	Net weight	175Kg
	Lifting range	480-630mm
	Net weight	43Kg

Designs and specification can be changed without prior notice for the purpose of improvement.

Ningbo Ming Sing Optical R&D Co.,Ltd

Tel:+86 574 87305541

E-mail: info@nbmingsing.com

Website: www.i-optik.com

Address: No.365 Middle Jingu Road(west), Ningbo, Zhejiang, China