



SHANGHAI PAX MEDICAL INSTRUMENT CO., LTD.  
Add: 3F B9 Building No.9 Lane 449, Nuijiang North Rd., Shanghai, China  
P.C.: 200333  
Tel: +86-21-52820699  
Fax: +86-21-52820763  
info@panalex.com.cn  
www.panalex.com.cn



panalex  
LED OPERATING LIGHT





- Using the medical class LED which is imported from the top lighting expert, German OSRAM Company.
- The light can satisfy different light spot requirements through its special optical software design.
- The unique spectral optical filter ensures the excellent color restoration.
- With over 50,000 hours effective working life, PANALEX can extremely save the hospital's maintenance cost.
- The energy-saving LED element with lowest heat emission features and cool light characteristic ensures the low surgical field temperature rise.

- PANALEX is easy for cleaning and can save more space through its streamline and ultrathin appearance design.
- The excellent lighting depth and shadowless effect makes the operations under different lighting depth more convenient.
- The detachable handle made of special material can be sterilized in 135 °C autoclave more than 3000 times.
- With 360° rotation and light arm, PANALEX ensures the light to be moved easily and freely.







PANALEX DVR

Digital Video Recorder

The system can control HD camera, indoor panoramic camera; The video information can be recorded, stored and broadcast to achieve remote live broadcasting and teaching. The storage capacity is 1000G.



Full HD Digital Camera System



1080p---Full HD Camera System

- The highest standard in medical field
- 1920\*1080 full HD resolution
- 16:9 digital image display
- Top level digital camera with best quality

HD Camera Specifications	
Image Sensor	1/2.8" Exmor CMOS
HD Digital Pixel(Resolution)	> 2 million pixels(1920×1080)
Aspect Ratio	16 : 9
Video Recording	1080p
Lens	10×optical zoom & 12×digital zoom Max optical zoom 120×
Aperture	F1.8-F3.4
Focal Length	f=3.8-38mm
Minimum Luminance	Colour 2.0 Lux, BlackWhite 0.1 Lux
White Balance	Indoor/fluorescent light/halogen light/LED light model, signal can be adjusted
Minimum Object Distance	800mm
Signal/ Noise Ratio	>50dB
Output Connection	1.HD SDI 2.HDMI



with emergency power



Technical Data	PANALEX
Central Illuminance (Lux)	160,000
Color Rendering Index(Ra)	>95
Color Temperature(K)	A:4500 B:3500 4000 4500 5000(Optional)
Focusable Light-field Size(mm)	150-300
Effective Bulb Service Life(hour)	>50,000
Diameter Lamp Housing(mm)	700
Lighting Depth L1+L2(mm)	>1000
Power(W)	80
MIS Endoscopy Lighting	Available

Options    ■ Remote control    ■ Touch-screem    ■ HD camera system