


Saturno-led

This is a surgical-type lamp suitable for minor surface operation surgeries, gynaecology and emergency room. The fact that the beams are close together (reflector size 195mm), means that they do not have to be focused. The lamp is very easy to move thanks to the lightness of the aluminium support structure.



SATPIN-LED



SATSONX2-LED



SATPAN-LED



SATSON-LED



Performances

Light intensity at 1 m distance (Ec)

Light head diameter

Color temperature

Color rendering index (CRI)

Diameter adjustment

d10 light field diameter where illuminance reached 10% of Ec

Light field diameter adjustable from-to

Depth of illumination IEC 60601-2-41 (L1+L2) at 60%

Depth of illumination IEC 60601-2-41 (L1+L2) at 20%

Total radiated energy Ee where the illuminance reaches max level

Ratio between radiated energy Ee and illuminance Ec

Average led life

Control of the illuminance

Electrical absorption

Saturno - led

50.000 lux

19,5 cm

2 selections: 4.000 - 4.500K

95 Ra

Fixed

260 mm

//

1100 mm

1800 mm

186 W/m²

3,63

> 60.000

20 - 100 %

19 W / 40 VA

Directive 93/42/EEC (main directive and further amendments) - Norm IEC 60601-2-41

ACCESSORIES



MAIN FEATURES

